













# MICROSTREAM™ CAPNOGRAPHY LINES

Microstream™ disposable CO<sub>2</sub> sampling lines provide continuous feedback on the breathing and ventilatory status of patients in emergency care. The lines help ensure accurate sampling for short term use, long term use, adult, pediatric and neonate care.



## Disposable CO<sub>2</sub> Sampling Lines

	Product	Part Number	Features	Connector	Compatibility
	Smart CapnoLine Plus with O <sub>2</sub> Connector, Adult	177653	Non-intubated, Oral-Nasal, with O <sub>2</sub> Delivery	 Orange: Microstream Connector for Short Term Use	 Philips MRx
	Smart CapnoLine Plus with O <sub>2</sub> Tubing, Adult	177268	Non-intubated, Oral-Nasal, with O <sub>2</sub> Delivery		
	Smart CapnoLine, Pediatric	177266	Non-intubated, Oral-Nasal		
	Smart CapnoLine with O <sub>2</sub> Tubing, Pediatric	177669	Non-intubated, Oral-Nasal, with O <sub>2</sub> Delivery		
	FilterLine Set, Adult/Pediatric	174620	Intubated, Filter		
	FilterLine Nasal Cannula, Infant	2722-47624	Non-Intubated, Nasal Sample		
	FilterLine H Set with Airway Adapter, Infant	176324	Intubated, Humidity Control For ETT 4.5 or smaller, 0.5 cc of deadspace	 Yellow: Microstream Connector for Long Term and Neonatal Use	 Lifepak 12 Physio-Control Lifepak 12 and 15 Philips MRX, MP2 Zoll X-Series
	FilterLine H Set Airway Adapter	M1921A	Intubated, Humidity Control		