

Bound Tree Medical Recertified AEDs and Monitor/Defibrillators

Monitoring & Defibrillation

AEDs and Monitors



Zoll M Series



LifePak 12



Philips MRx



PROFESSIONAL AEDs

For more information contact Bound Tree Medical at 800.533.0523 or search our online inventory of recertified equipment at www.BoundTree.com.

In today's economic climate, more and more decision makers are considering the purchase of used equipment in lieu of new. To help ensure that our customers receive the best equipment possible, Bound Tree Medical offers a recertified option for AEDs and Monitor/Defibrillators.

There are many companies that sell refurbished equipment. But what does "refurbished" actually mean? Is the device simply wiped down and turned on to see if it is functional? Or is it put through the same rigorous testing that it goes through at the manufacturer, by factory-certified technicians, to ensure it meets the exact specifications that were required of it when it came off the assembly line?

State-of-the-art Testing

We only use current, state-of-the-art testing equipment from Fluke and BCI to scrutinize how each parameter of your device performs. These devices are regularly calibrated and maintained according to the factory recommendations. Your tested device will include a calibration sticker attached to the device and a results sheet for your records at your request.

Detail, Detail, Detail

Your device is painstakingly detailed to ensure your device is as close to out-of-the-box new in appearance as possible. We then put only new electrodes and batteries on every device (something not all suppliers do).

We Have You Covered

We carry full liability coverage that extends to the equipment we sell (although we have never needed to file a claim). We take the recertification process seriously and you can count on us to treat your device as if our own lives depended on it.

Not all used equipment is alike. Selling pre-owned medical equipment is a serious business and deciding on which company to purchase from is a vital decision. Demand the highest quality recertified equipment backed by Bound Tree's unrivaled service. If peace of mind and premium recertified equipment are priorities of yours, Bound Tree Medical continues to be your partner in EMS.

BoundTree
medical

800-533-0523 | www.boundtree.com

Our 11-point checklist for AEDs, and 22-point checklist for monitor/defibrillators, are made by trained technicians in order to adhere to stringent manufacturer guidelines and specifications.

11-Point AED Checklist

1. AED will be upgraded to current AHA guidelines (if upgrade available from manufacturer)
2. Each unit is connected to a new Fluke Impulse 7000 defibrillator analyzer to simulate a shockable rhythm, and 3 shocks are delivered to assure the joule output is within FDA guidelines
3. All device and performance information is recorded and stored in both electronic and paper form
4. Visual inspection to ensure AED is in acceptable cosmetic condition
5. Device is thoroughly and meticulously cleaned (think Q-tips and toothpicks)
6. AED is accessorized and made patient-ready with electrodes and batteries installed
7. AED is turned on to initiate self-test that assures it is ready to be used in a rescue
8. Battery and Pad expiration dates are recorded to ensure you are notified 60 days before your accessories require replacement (end user ultimately responsible for replacement)
9. Packing slip is hand 'checked' to ensure all items are included
10. Order is professionally packed in accordance to Fedex/UPS regulations
11. ED tracing is submitted to the manufacturer in accordance with FDA regulations

22-Point Monitor/Defibrillator Checklist

1. Unit is visually inspected to ensure device is cosmetically acceptable
2. Software version is recorded
3. Software is upgraded to latest version (if available from manufacturer)
4. Time and Date are set
5. Unit and customer information is entered into Ansur Software and stored in both electronic and paper form
6. Device is tested on equipment regularly calibrated by the manufacturer
7. Joule output is tested and recorded at intervals of: 10, 20, 50, 100, 200, 300, 360 (if applicable)
8. Joule output is recorded and displayed on a sticker placed on the device
9. Calibration sticker notes "Next Inspection" date according to manufacturer guidelines
10. Pacing is tested (if applicable)
11. Spo2 is tested (if applicable)
12. Non-invasive Blood Pressure is tested (if applicable)
13. Capnography is tested (if applicable)
14. Printer is tested, and test strip is shipped with the unit to show functionality
15. All testing results stored on server in case ever required and test summary is printed and filed
16. Unit is fully accessorized and assembled patient-ready
17. Packing slip is printed and hand checked it insure all items are included
18. Device is thoroughly and meticulously cleaned (think Q-tips and toothpicks)
19. Device is professionally packed using approved materials and insured
20. Tracking # is recorded on SO/Invoice
21. Battery and Pad expiration dates are recorded to ensure you are notified 60 days before your accessories require replacement (end user ultimately responsible for replacement)
22. Device information and destination are submitted to the manufacturer in accordance with FDA regulations