


MedSource

ClearSafe™ Safety I.V. Catheter

I.V., Drug Delivery &
Blood Collection

I.V. Catheters

11.2015



ClearSafe™ retracted

Ordering Information

All are available in 1/EA, 50/BX, or 200/CS

ClearSafe™ Safety I.V. Catheter

1612-84110	14 ga x 1 1/4"
1612-84120	16 ga x 1 1/4"
1612-84130	18 ga x 1 1/4"
1612-84140	20 ga x 1 1/4"
1612-84145	20 ga x 1"
1612-84150	22 ga x 1"
1612-84160	24 ga x 3/4"

ClearSafe™ Comfort I.V. Safety Catheter

1612-84210	14 ga x 1 1/4"
1612-84220	16 ga x 1 1/4"
1612-84230	18 ga x 1 1/4"
1612-84240	20 ga x 1 1/4"
1612-84245	20 ga x 1"
1612-84250	22 ga x 1"
1612-84260	24 ga x 3/4"

ClearSafe™ Safety I.V. Catheter

The cost-effective MedSource ClearSafe™ I.V. catheter is designed to allow caregivers easy and safe I.V. access. The specially designed beveled needle provides smooth, painless insertion. The clicking lock mechanism ensures safety from cross-contamination.

- Fully encapsulating needle shield with clicking lock safety mechanism
- Ergonomically designed thumb grip
- Gently tapered catheter tip
- Large-volume chamber
- Designed to help healthcare workers protect themselves from exposure to blood borne pathogens via needlestick injuries
- Provides smooth transition from needle to catheter
- Allows for a variety of user techniques to aid insertion
- Provides clear visual conformation through flashback

ClearSafe Comfort™ Safety I.V. Catheter

The ClearSafe Comfort contains all of the features and benefits as its counterpart, the ClearSafe, but with the added benefit of PUR (polyurethane) catheter material. This material provides a softer, more comfortable patient experience, and a longer in-dwelling performance for the catheter.

For more information or a product demonstration please contact your local Account Manager or visit www.boundtree.com/contactus.aspx

BoundTree
medical

800.533.0523 | www.boundtree.com