

Onyx[®] Vantage 9590 Finger Pulse Oximeter



The Onyx Vantage 9590 finger pulse oximeter with PureSat[®] technology can quickly and accurately capture SpO₂ and pulse rate measurements, even on patients where low perfusion is a challenge. The Onyx comes with lanyard, instructions on CD and two AAA batteries. Optional carrying case available.

- Oxygen saturation range — 0% to 100%
- Pulse rate range — 18 to 300 pulses/minute
- Comfortable to carry at just two ounces
- Easy battery access
- Lead-free and latex-free construction
- Made in the USA

One Year Free Insurance program for the Onyx Vantage 9590 if it is lost or stolen within one year from the date of purchase, Nonin will provide a one-time free replacement unit. To enroll in insurance, register your serial number at: www.nonin.com/insurance

Ordering Information

Onyx Vantage 9590
2711-89590

Justice Mark II, Hard Case for the Onyx
173802

Product Features

- Accurate — Wide range of patient populations and settings, including low perfusion and motion
- Durable — Resistant to dropping, water spills/ingress and reads in challenging environments where it can go from hot to cold and cold to hot temperatures
- Efficient — Up to 6,000 spot checks or 36 hours of continual operation on two AAA batteries
- Versatile — Works on pediatric to large adult patients, with claims for use on all fingers, thumbs, toes & dark skin tones
- Available in four colors — Black, blue, red and purple
- Warranty — 4-year warranty
- Wide range of finger thicknesses — From 8 mm to 25.4 mm (0.3" to 1.0")
- Clinician-facing, easy-to-see, bright LED — SpO₂ and pulse rate can be read at an angle in bright or low-light conditions
- Three-color pulse quality indicator LED — Provides a quick assessment of patient pulse quality

For more information, please contact your Account Manager or call Customer Service at 800-533-0523.

