

FEATURED TECHNOLOGY

ADAPTABLE C-CURVE™ PRINCIPLE

SAM® Soft Shell Splint leverages the revolutionary C-Curve™ Principle to mold its foam and aluminum layers into structural arcs, strong enough for pre- or post-cast care, while remaining lighter and more transportable than traditional splints that rely on heavy, rigid materials. Bend it back into its original form and it's ready to be remolded.





SAM Soft Shell

ENGINEERED FOR VERSATILE SUPPORT

MAXIMUM-COMFORT SURFACE

Absorbent-brushed terry cloth fabric covers skin-contact surface, providing enhanced comfort.

2 MOLDABLE ALUMINUM

Flat O-temper aluminum inside the splint's soft outer layer molds easily, allowing for a wide spectrum of applications.

3 RADIOLUCENT MATERIAL

Medical professionals are able to take X-rays and CT-scans without removing the splint from the fractured or injured limb.

SIZE THAT FITS

Full product line covers most sizes needed, including XL and sugar tong.


RESISTANCE CONTROL EDGE

Firm and extra-firm ends allow different levels of resistance to motion, according to needs. (Not available on the Sugar Tongs.)

FIRM

OR

EXTRA FIRM



MODELS







XL 5.5 in wide





SUGAR TONGPediatric, Adult, Adult XL

PRODUCT SPECIFICATIONS

MODELS AVAILABLE*

EC REP

9" (SS540-EN) 9" XL (SS541-EN) 12" (SS542-EN) 12" XL (SS543-EN) 15" (SS544-EN) 15" XL (SS545-EN)

X-RAY RADIOLUCENT



Emergo Europe, Prinsessegracht 20, 2514 AP, The Hague, The Netherlands SAM Medical Products®
12200 SW Tualatin Rd., Ste 200
Tualatin, OR 97062 USA

SUGAR TONG MODELS

Pedatric (ST561-EN) Adult (ST560-EN) Adult XL (ST562-EN)

GLOBAL SALES

U.S. Toll Free: 800.818.4726 U.S. Local: 503.639.5474 U.S. Fax: 503.639.5425 **sammedical.com**

* Standard width 4.25 in wide, XL 5.5" wide.

SAM® is a registered trademark in the United States and other countries.

©2019 SAM Medical Products. All rights reserved.